

NUTRICAIR™

BABY BOTTLES

β sterilized



*Pediatric enteral nutrition
without any compromise*



In the beginning, there was a patient...



The guarantee of a sterile product secured against infectious and toxic risks without any compromise

- » Feeding bottles, containers and teats suitable for all newborns, without risk of exposure to toxic gas residues₍₁₎₍₂₎₍₃₎
- » Assurance of sterility: 10⁻⁶ ufc/ml in accordance with EMA statement₍₄₎
- » CE marked and compliant with international standards ISO 10993-1, ISO 11137-1 and EU regulation n°10/2011
- » Environment friendly sterilization process₍₂₎₍₅₎

Cair LGL commitment

Because our products are intended for prematured babies, newborns, infants with pathologies, in a medical environment,

Because we do not want these children to be exposed to potential pathogens,

Because for us, providing a non-sterile or bacteriologically clean bottle to a baby in hospital is unsatisfactory,

» we have made the choice, for Nutricair™ bottle range, to offer only

sterile products

Because having to make a benefit/risk trade-off between the sterilisation method used and the infectious and toxic risks is not, for us, a satisfactory long-term solution,

And because our duty is to do everything possible to respect the recommendations of the DGS₍₆₎

» we have made the choice, for Nutricair™ bottle range,

sterilisation by beta radiation

Prepare, preserve and feed: the bottles

» Easy to handle

» Easy to screw and unscrew

Bottles with teat

- Preparation, administration
 - » Teats : 3 speeds for prematured babies, 1 speed, 3 speeds
 - » Bottles : 60, 130 and 240 ml



Bottles with cap

- Collection, preparation, storage
 - » Flexible and robust polypropylene graduated bottles
 - » Waterproof, suitable for pasteurization and freezing



Sucking: teats

» 3 teats adapted to breast milk

» for 3 infant sucking modes

Prematured babies

3 speeds



Full-term infant

1 speed

3 speeds



- Designed for comfortable sucking

- » Soft elastomer teats
- » Regulation system of the air inlet in the bottle

- Immediately recognizable

- » A color code for each type of teat

- Respect of the hygiene

- » Presented with teat cap: graduated milk cup marked from 10 to 60 ml and clipped on the ring

Milk cup

- Graduated from 10 to 60 ml

